

Unitol®

the separator

Isopraim®

one separator multiple applications



PERIDENT DENTAL PRODUCTS s.r.l.

Via del Fornaccio, 46/F - località Vallina - 50012 Bagno a Ripoli
FIRENZE (ITALIA) Tel. (+39) 055 69 65 40



PERIDENT DENTAL PRODUCTS

www.perident.it

info@perident.it

only one separator for multiple applications

Isopraim is a multipurpose silicone-based separator with high degree of separation.

This products boasts 7 specific applications, granting top results to selective dental technicians.

To separate :

- gypsum - gypsum
- gypsum - wax
- gypsum - silicone
- wax - wax
- wax - metal
- silicone - silicone
- silicone - polyether



Separation of flasks and creation of pindex models



Impressions pouring with gingival reproduction



Modelling copings of preformed abutement



REF 107 - Isopraim 280 ml bottle

benefits :

- optimum spreadability
- creates no thickness
- dries quickly
- leaves no trace
- excellent detachability
- Alcohol and solvents free

the separator for plaster - resin

Unifol is an insulating containing with alginate base and emulsified oils with exceptional separating properties.

Used on both hot or cold moulds, it creates an effective barrier between plaster and resin, avoiding a bleaching effect.

To separate :

- gypsum - resin



REF 101 - unifol 450 ml bottle
REF 103 - unifol triple pack (3 bottles x 450 ml)

Benefits :

- reduces labor when deflasking
- reduces finishing and polishing time
- leaves no residual film and does not create thickness
- leaves surfaces neat and clean



Flask curing and injection curing techniques



Cold resins and composites



Light-curing materials