

Unitrol®

The separator



PERIDENT DENTAL PRODUCTS s.r.l.
Via del Fornaccio, 46/F - località Vallina - 50012 Bagno a Ripoli
FIRENZE (ITALIA) Tel. (+39) 055 69 65 40



PERIDENT
dental products

www.perident.it

info@perident.it


The separator for plaster - resin

unifol® is an insulating containing a base alginate and emulsified oils, producing exceptional separating properties.

Used on both hot or cold moulds, it creates an effective barrier between plaster and resin, avoiding a bleaching effect.

Benefits:

- Reduces labor when deflasking
- Reduces finishing and polishing time
- Leaves no residual film to take up space
- Restorations exhibits smoother and cleaner surfaces



Flask curing and injection curing techniques



Cold resins and composites



Light-curing materials



REF. 101 - **unifol®** 450 ml bottle

REF. 103 - **unifol®** triple pack (3 bottles x 450 ml)