

# Unitrol<sup>®</sup>

*The separator*



PERIDENT DENTAL PRODUCTS s.r.l.  
Via del Fornaccio, 46/F - località Vallina - 50012 Bagno a Ripoli  
FIRENZE (ITALIA) Tel. (+39) 055 69 65 40



**PERIDENT**  
dental products

[www.perident.it](http://www.perident.it)

[info@perident.it](mailto:info@perident.it)

# *The separator for plaster - resin*

**unifol®** is an insulating containing a base alginate and emulsified oils, producing exceptional separating properties.

Used on both hot or cold moulds, it creates an effective barrier between plaster and resin, avoiding a bleaching effect.


## *Benefits:*

- Reduces labor when deflasking
- Reduces finishing and polishing time
- Leaves no residual film to take up space
- Restorations exhibits smoother and cleaner surfaces



REF 101 - **unifol®** 450 ml bottle

REF 103 - **unifol®** triple pack (3 bottles x 450 ml)



Flask curing and injection curing techniques



Cold resins and composites



Light-curing materials